

Made in Germany predominantly for Law Enforcement applications, the Dräger DrugCheck® 3000 is now available in Australia for Medical, WorkPlace or Personal Drug Testing - this is a Police grade saliva drug test.

Dräger DrugCheck® 3000

Obtain samples easily and safely

The Dräger DrugCheck 3000 consists of two components: an absorbent collector for obtaining a saliva sample, and a test cassette for the analysis. The test cassette contains the buffer liquid and a window with two test strips, which displays the control and test lines.

The drug test itself is performed in three easy steps: swab the saliva sample, shake the test kit, wait for a brief incubation period, and then start the test. A colour indicator on the swab disappears as soon as it has absorbed enough oral fluid for a test.

As soon as the control lines appear in the window, you can read the results. If the test result is negative, a line will appear alongside the respective substance class (drug). This means that none of the target substance was detected in the sample. If a line fails to appear next to one of the substance classes, then the result for this substance is positive.

Fast and sensitive drug analysis for up to 6 substance classes

Check individuals for up to six substance classes simultaneously with Dräger DrugCheck 3000: **cocaine, opiates, amphetamine, methamphetamines/designer drugs (e.g. Ecstasy, MDMA), benzodiazepines (e.g. in prescription medicines) and cannabis (THC)**. Of all the substance classes listed, cannabis is the drug consumed most frequently and also the most difficult of all compounds to identify (THC = delta-9-tetrahydrocannabinol).



Now detects up to 6 substance classes

- » Detects up to 6 substance classes
- » High sensitivity (THC down to 15ng/ml)
- » High accuracy
- » Quick sample collection <15 seconds
- » Fast analysis time of 3 minutes
- » Hygienic and intuitive handling
- » Operating temperature 5 - 40°C
- » Manufactured in Germany

This is why the Dräger DrugCheck 3000 was optimised to detect THC, and now offers two measurement options: depending on the chosen waiting time (pre-incubation) of the test there is a fast option with a higher or sensitive option with a low THC cut-off. The entire test only takes 3 – 5 minutes.

Unambiguous on-the-spot drug detection


Its compact design makes the DrugCheck® 3000 drug test easy to transport. It can be made ready for use quickly and easily, and is suitable for roadside drug checks and drug testing in safety-relevant workplaces.

The test kit has no electrical parts, which makes the drug tester safe to use in hazardous areas and even as point-of-care-testing (POCT).

The disposable DrugCheck 3000 test kit cannot be manipulated and is hygienic to use.



Features



Scan to watch VIDEO!



Scan the QR Code using an iOS / Android App to watch a video on the *DrugCheck® 3000!*

Accessories

The Dräger DCD 5000 allows for an effective method of collecting and storing saliva samples that are to be sent to the laboratory for further analysis – without the sample being falsified along the way.



Technical Data

| | |
|--------------------------------|---|
| Dimensions (W x H x D) | 32 x 111 x 57 mm 1.26 x 4.37 x 2.24 inches |
| Weight | < 30 g (0,07 lb) |
| Operating range | Operation: 5 °C to 40 °C (41 °F to 104 °F) at 5 % to 95 % relative humidity |
| Storage / transportation | 4 °C to 25 °C (39 °F to 77 °F) |
| Time of measurement | Fast measurement: analysis < 3 minutes Sensitive measurement: analysis < 5 minutes |
| Selecting the measurement mode | Depending on the selected pre-incubation time and the desired THC cut-off fast or more sensitive reading. |
| Licenses / standards | Licensed as a medical product within the EU in accordance with Directive 98/79/EC on in vitro diagnostic medical devices (IVD). Outside of the European Union for non-medical use (Non-IVD) or for forensic use (only USA). |
| Cut off Level | THC 15 or 40 ng/ml, COC 20 ng/ml, OPI 20 ng/ml, AMP 50 ng/ml, MET 50 ng/ml, BZ 15 ng/ml. |

